




| | | |
|---|--|--|
| theranos | Anti-Streptolysin O (ASO) Validation Report | Document Number: CL-RPT-145002 |
| | Validation Document | Revision: A Effective Date: 1/17/2014 |
| Antistreptolysin O (ASO) Assay Validation Report | | |

Author(s):


| | |
|--|----------------------------|
| Signature:  | Date: 2/19/2014 |
| Name: Adam Rosendorff, MD | Title: Laboratory Director |

Reviewer(s)

| | |
|--|--------------------------------|
| Signature:  | Date: 3/19/14 |
| Name: Sharada Sivaraman, PhD | Title: Immunoassay Team Leader |

| | |
|---|-----------------------|
| Signature:  | Date: 3/18/2014 |
| Name: Daniel Young, Ph.D. | Title: Vice President |

Approver(s):

| | |
|--|----------------------------|
| Signature:  | Date: 2/19/2014 |
| Name: Adam Rosendorff, M.D | Title: Laboratory Director |

 2/19/15
Sunil S. Dhawan M.D.

| | | |
|---|--|---|
| theranos | Anti-Streptolysin O (ASO) Validation Report | Document Number: CL-RPT-145002 Revision: A |
| | Validation Document | Effective Date: 1/17/2014 |
| Antistreptolysin O (ASO) Assay Validation Report | | |

TABLE OF CONTENTS

| | | |
|---|--|------|
| 1 | ASSAY BACKGROUND | 3 |
| 2 | REGULATION AND GUIDANCE | 3 |
| 3 | PRINCIPLE OF THE PROCEDURE | 3 |
| 4 | CALIBRATION | 4 |
| 5 | QUALITY CONTROL | 5 |
| 6 | PRECISION | 5-8 |
| 7 | METHOD COMPARISON: THERANOS VERSUS PREDICATE | 8-10 |
| 8 | REFERENCE RANGE VERIFICATION | 10 |
| 9 | REFERENCES | 11 |

| | | |
|---|--|--------------------------------|
| theranos | Anti-Streptolysin O (ASO) Validation Report | Document Number: CL-RPT-145002 |
| | Validation Document | Revision: A |
| Antistreptolysin O (ASO) Assay Validation Report | | Effective Date: 1/17/2014 |

1 ASSAY BACKGROUND

During acute infection with group-A hemolytic Streptococci, antibodies to bacterial exotoxins are produced, that react with streptolysin O coated latex beads. ASO reactivity is also seen in rheumatic fever (valvular heart disease), and streptococcal glomerulonephritis. ASO can be used as an adjunct to the rapid strep test in the acute care setting.

2 REGULATION AND GUIDANCE

The qualification/validation of the ELISA assays on the Theranos device will be in accordance with C.F.R. Ch IV, § 493.1253 "Standard: Establishment and verification of performance specifications" and outlined in CLSI guideline C28A3.

3 PRINCIPLE OF THE PROCEDURE

In those infections promoted by acute streptococcal infection, antibodies to the exotoxin of streptococcus are usually produced. By reacting suspended uniform latex particles coated with streptolysin-O with serum containing antibodies, an increase in turbidity occurs. By comparing with a standard, a quantitative value for the concentration of anti-streptolysin-O present in serum is obtained.

Plasma samples were diluted 1:3,125 fold in saline prior to analysis.

| | | |
|---|--|--|
| theranos | Anti-Streptolysin O (ASO) Validation Report | Document Number: CL-RPT-145002 |
| | Validation Document | Revision: A Effective Date: 1/17/2014 |
| Antistreptolysin O (ASO) Assay Validation Report | | |

4 CALIBRATION

- 4.1 In 42 CFR Part 493.1255, it is required to perform calibration procedures with at least the frequency recommended by the manufacturer, or using criteria specified by the laboratory, or when calibration verification fails to meet acceptable limits.
- 4.1.1 The term "calibration verification," as used in CLIA, includes:
- 4.1.1.1 Confirming that a calibration meets the method manufacturer's specifications
- 4.1.1.2 Verifying that the calibration is suitable for the entire measuring interval (or "reportable range," which is the CLIA term)
- 4.2 Calibrators were diluted 1:3.125 and verified on the ADVIA system
- 4.2.1 This dilution factor is within the acceptable limits of the ADVIA internal calibration test.
- 4.3 For the purposes of this Validation Plan, calibration was carried out with every new lot of reagents.
- 4.3.1 Each level was tested in replicates of 3 and the average was used to create a standard curve for testing.
- 4.3.2 The calibration was verified using quality controls.

| | | |
|---|--|--|
| theranos | Anti-Streptolysin O (ASO) Validation Report | Document Number: CL-RPT-145002 |
| | Validation Document | Revision: A Effective Date: 1/17/2014 |
| Antistreptolysin O (ASO) Assay Validation Report | | |

5 QUALITY CONTROL

5.1 Two to four level quality control samples, as appropriate to the assay, were analyzed with each calibration and before each test during the validation.

5.1.1 High = 85 IU/mL

5.1.2 Mid = 156 IU/mL


5.1.3 Low = 223 IU/mL

5.2 The QC levels are not included when generating the calibration curve.

6 PRECISION

6.1 Precision was evaluated according to CLSI standard EP5-A2, Evaluation of Precision Performance of Quantitative Measurement Methods.

6.1.1 A total of 20 runs were performed over 10 days with 2 runs per day and 2 replicates per run for a total of 40 data points. The following tables indicate the between-run, between-day and within-laboratory precision at 3 levels (High = 223 IU/mL, Mid = 156 IU/mL, Low = 85 IU/mL). The following data describes the results obtained:

| | | |
|---|--|---|
|  | Anti-Streptolysin O (ASO) Validation Report | Document Number: CL-RPT-145002 Revision: A |
| | Validation Document | Effective Date: 1/17/2014 |
| Antistreptolysin O (ASO) Assay Validation Report | | |

Precision

CLSI guideline EP05-A2 section 10.8

Level = L1

| | |
|------------------------|----|
| Number of observations | 40 |
| Number of runs | 20 |
| Number of days | 10 |
| Runs per day | 2 |
| Replicates per run | 2 |

CLSI guideline EP05-A2 section 10.4 recommends a minimum of 40 runs, with 2 replicates per run.

| | | | | |
|-------------------|-------|--------------|------|--------------------|
| Mean | 81.99 | | | |
| | SD | 95% CI | CV | Allowable Total SD |
| Repeatability | 1.89 | 1.45 to 2.73 | 2.3% | - |
| Between-run | 0.00 | | 0.0% | - |
| Between-day | 1.66 | | 2.0% | - |
| Within-laboratory | 2.52 | 1.96 to 3.53 | 3.1% | 20.50 |

Imprecision is less than allowable total imprecision: 25%.

Level = L2

| | |
|------------------------|----|
| Number of observations | 40 |
| Number of runs | 20 |
| Number of days | 10 |
| Runs per day | 2 |
| Replicates per run | 2 |

CLSI guideline EP05-A2 section 10.4 recommends a minimum of 40 runs, with 2 replicates per run.

| | | | | |
|-------------------|--------|--------------|------|--------------------|
| Mean | 153.11 | | | |
| | SD | 95% CI | CV | Allowable Total SD |
| Repeatability | 2.12 | 1.62 to 3.06 | 1.4% | - |
| Between-run | 2.81 | | 1.8% | - |
| Between-day | 1.41 | | 0.9% | - |
| Within-laboratory | 3.79 | 2.97 to 5.24 | 2.5% | 38.28 |

Imprecision is less than allowable total imprecision: 25%.

| | | |
|---|--|--|
| theranos | Anti-Streptolysin O (ASO) Validation Report | Document Number: CL-RPT-145002 |
| | Validation Document | Revision: A Effective Date: 1/17/2014 |
| Antistreptolysin O (ASO) Assay Validation Report | | |

Level = L3

| | |
|------------------------|----|
| Number of observations | 40 |
| Number of runs | 20 |
| Number of days | 10 |
| Runs per day | 2 |
| Replicates per run | 2 |

CLSI guideline EP05-A2 section 10.4 recommends a minimum of 40 runs, with 2 replicates per run.

| | | | | |
|-------------------|--------|--------------|------|-----------------------|
| Mean | 223.15 | | | |
| | SD | 95% CI | CV | Allowable Total SD |
| Repeatability | 2.80 | 2.14 to 4.04 | 1.3% | - |
| Between-run | 3.33 | | 1.5% | - |
| Between-day | 3.74 | | 1.7% | - |
| Within-laboratory | 5.73 | 4.35 to 8.40 | 2.6% | 55.79 |

Imprecision is less than allowable total imprecision: 25%.

Theranos

| | | |
|---|--|---|
| theranos | Anti-Streptolysin O (ASO) Validation Report | Document Number: CL-RPT-145002 Revision: A |
| | Validation Document | Effective Date: 1/17/2014 |
| Antistreptolysin O (ASO) Assay Validation Report | | |

6.1.2 The mean recovery was as follows:

| <u>Control#</u> | <u>Assigned (mg/dL)</u> | <u>Theranos (mg/dL)</u> | <u>% Recovery</u> |
|-----------------|-------------------------|-------------------------|-------------------|
| 1 | 85 | 82 | 96.4% |
| 2 | 156 | 153 | 98% |
| 3 | 223 | 223 | 100% |

6.2 Acceptance criteria:

Total allowable error (TAE %) of 25%, was selected as the acceptance criteria for this assay following CLIA proficiency guidelines as printed in the Federal Register February 28, 1992;57(40):7002-186. SD as well as total error guidelines as stated in the American Proficiency Institute Peer Data for 2013 CHEMISTRY / IMMUNOLOGY / IMMUNOHEM -1ST EVENT. Allowable bias was calculated as the residual error budget after precision values (CV %) were subtracted from TAE (%). Values in brackets indicate the closest/corresponding API levels obtained from 2013 proficiency data.

Table II

| | Level 1 (82) | Level 2 (153) | Level 3 (223) |
|----------------------|---------------------|----------------------|----------------------|
| TAE% | 25 | 25 | 25 |
| CV (%) | 3.1 | 2.5 | 2.6 |
| Allowable Bias (%) | 21.9 | 22.5 | 22.4 |
| Bias (%) (Table III) | 9 | 9 | 9 |

7 METHOD COMPARISON, THERANOS VERSUS PREDICATE

7.1 Twenty (20) venous samples were run using the predicate Siemens protocol without dilution, and in parallel on the Theranos assay with pre-dilution. Results were plotted in a scatter diagram, and a simple linear regression was performed (Figure I). Raw data as well as the scatter-plot summarizing the results are shown in Table III.

7.2 Mean bias comparing methods was calculated as follows: %Bias=[(Theranos-Siemens)/Siemens]*100 and results are shown in the column labelled "% difference" and indicated in Section 6.2.

7.3 Mean bias is less than allowable bias therefore, the acceptance criteria PASS.

| | | |
|---|--|--|
| theranos | Anti-Streptolysin O (ASO) Validation Report | Document Number: CL-RPT-145002 |
| | Validation Document | Revision: A Effective Date: 1/17/2014 |
| Antistreptolysin O (ASO) Assay Validation Report | | |




Figure 1: Bias estimation, Theranos versus Siemens assay

| Sample # | Siemens | EDTA | Theranos | EDTA | % difference | T-corr | %difference |
|----------|---------|-------|----------|-------|--------------|--------|-------------|
| 1 | ASO_2 | 398.0 | T-ASO | 398.0 | 0% | 396.0 | 0.00% |
| 2 | ASO_2 | 244.4 | T-ASO | 235.5 | 4% | 234.6 | 3.64% |
| 3 | ASO_2 | 166.0 | T-ASO | 170.2 | -2% | 169.7 | -2.53% |
| 4 | ASO_2 | 151.9 | T-ASO | 143.5 | 6% | 143.1 | 5.53% |
| 5 | ASO_2 | 151.9 | T-ASO | 146.4 | 4% | 146.0 | 3.62% |
| 6 | ASO_2 | 149.3 | T-ASO | 142.6 | 5% | 142.2 | 4.49% |
| 7 | ASO_2 | 117.3 | T-ASO | 112.1 | 5% | 111.9 | 4.43% |
| 8 | ASO_2 | 117.0 | T-ASO | 113.8 | 3% | 113.6 | 2.74% |
| 9 | ASO_2 | 116.6 | T-ASO | 115.6 | 1% | 115.4 | 0.86% |
| 10 | ASO_2 | 99.3 | T-ASO | 90.5 | 9% | 90.5 | 8.86% |
| 11 | ASO_2 | 84.3 | T-ASO | 76.5 | 10% | 76.5 | 9.25% |
| 12 | ASO_2 | 71.9 | T-ASO | 64.0 | 12% | 64.1 | 10.99% |
| 13 | ASO_2 | 62.7 | T-ASO | 58.9 | 6% | 59.1 | 6.06% |
| 14 | ASO_2 | 47.9 | T-ASO | 44.6 | 7% | 44.8 | 6.89% |
| 15 | ASO_2 | 41.1 | T-ASO | 33.0 | 22% | 33.3 | 19.71% |
| 16 | ASO_2 | 40.4 | T-ASO | 37.0 | 9% | 37.3 | 8.42% |
| 17 | ASO_2 | 35.4 | T-ASO | 30.7 | 14% | 31.0 | 13.28% |

| | | |
|---|--|--------------------------------|
| theranos | Anti-Streptolysin O (ASO) Validation Report | Document Number: CL-RPT-145002 |
| | Validation Document | Revision: A |
| Antistreptolysin O (ASO) Assay Validation Report | | Effective Date: 1/17/2014 |

| | | | | | | | |
|---------|-------|-------|-------|-------|------|-------|--------|
| 18 | ASO_2 | 28.8 | T-ASO | 23.0 | 22% | 23.4 | 20.14% |
| 19 | ASO_2 | 18.3 | T-ASO | 16.4 | 11% | 16.8 | 10.38% |
| 20 | ASO_2 | 11.1 | T-ASO | 7.9 | 34% | 8.4 | 28.83% |
| Average | | 107.7 | | 103.0 | 9.0% | 102.9 | 8.28% |

Table III Bias estimation, Theranos versus Siemens assay

8 CTN REFERENCE RANGE VERIFICATION

- 8.1 20 unique capillary tube and nanotainer (CTN) samples were collected from healthy donors and assayed in singlicate using the Theranos methods, as shown in Table III. Resulting values were corrected to match more closely with the predicate using the regression equation as follows: Corrected value=(CTN value +5.3516)/1.0063. The corrected values are shown in table III, under T-corr.
- 8.2 19/20 (95%) of corrected CTN values fell within the predicate reference range (<240 IU/mL), therefore the reference range is verified (CLSI C28-A3).

| | | |
|---|--|--------------------------------|
| theranos | Anti-Streptolysin O (ASO) Validation Report | Document Number: CL-RPT-145002 |
| | Validation Document | Revision: A |
| Antistreptolysin O (ASO) Assay Validation Report | | Effective Date: 1/17/2014 |

9 REFERENCES

- 9.1 Code of Federal Regulations, Title 42, Chapter IV, Subchapter G, Part 493, Subpart K, Sections 493.1217, 493.1253, and 493.1255.
- 9.2 DeSilva B, Smith W, Weiner R, et al. Recommendations for the bioanalytical method validation of ligand-binding assays to support pharmacokinetic assessments of macromolecules. *Pharmaceutical Res.* 2003; 20:1885-1900.
- 9.3 Guidance for Industry: bioanalytical method validation. U.S. Department of Health and Human Services, Food and Drug Administration, Center for Drug Evaluation and Research, 2001.
- 9.4 R (version 2.13.1). The R Foundation for Statistical Computing, 07/08/2011.
- 9.5 StatisPro (version 1.13.00). Clinical and Laboratory and Standards Institute, Wayne, PA. 07/14/2011.
- 9.6 Dexter-Immunoassay (version 1.0), Theranos, Inc., 2009.
- 9.7 EP10-A3, Preliminary Evaluation of Quantitative Clinical Laboratory Measurement Procedures; Approved Guideline—Third Edition, 2006, Clinical and Laboratory Standards Institute, Wayne, PA.
- 9.8 EP15-A2, User Verification of Performance for Precision and Trueness; Approved Guideline—Second Edition. 2005, Clinical and Laboratory Standards Institute, Wayne, PA.
- 9.9 EP09-A2-IR, Method Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Second Edition (Interim Revision), 2010, Clinical and Laboratory Standards Institute, Wayne, PA.
- 9.10 EP05-A2, Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline—Second Edition, 2004, Clinical and Laboratory Standards Institute, Wayne, PA.
- 9.11 EP06-A, Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline, 2003, Clinical and Laboratory Standards Institute, Wayne, PA.
- 9.12 EP7-A2, Interference Testing in Clinical Chemistry; Approved Guideline—Second Edition, 2004, Clinical and Laboratory Standards Institute, Wayne, PA.

| | | |
|---|--|---|
| theranos | Anti-Streptolysin O (ASO) Validation Report | Document Number: CL-RPT-145002 Revision: A |
| | Validation Document | Effective Date: 1/17/2014 |
| Antistreptolysin O (ASO) Assay Validation Report | | |

Theranos Internal
Only