



theranos	LDT Validation Report	Theranos Phosphorus Assay	Rev:
		CL RPT-14057	1
Description	Validation Report for Modified Siemens Assay of Inorganic Phosphorus in Lithium Heparin Plasma		
Originator: Curtis Schneider		Date: 10/15/2013	

Validation of Modified Siemens Inorganic Phosphorus Assay


Author(s):


Signature: 	Date: 12/12/13
Name: Paul Patel, Ph.D.	Title: Team Lead, General Chemistry

Reviewer(s):

Signature:	Date:
Name:	Title:
Signature: 	Date: 12/15/2013
Name: Daniel Young, Ph.D.	Title: Vice President

Approver(s):

Signature: 	Date: 12/12/2013
Name: Adam Rosendorff, M.D.	Title: Laboratory Director

 9/19/15
Sunil S. Dhawan M.D.

Theranos Confidential
Page 1 of 15

Any retransmissions, reproductions, dissemination or other use of these materials by persons or entities other than the intended recipient is prohibited. This document supersedes all earlier or previous documents unless approved in writing.

TMP-00009 Rev. A, Released 08/01/13

theranos	LDT Validation Report	Theranos Phosphorus Assay	Rev:
		CL RPT-14057	1
Description	Validation Report for Modified Siemens Assay of Inorganic Phosphorus in Lithium Heparin Plasma		
Originator: Curtis Schneider		Date: 10/15/2013	

Inorganic Phosphorus Plasma Assay

- I. **Overview**
- II. **Method Principle**
- III. **Definitions and Abbreviations**
- IV. **Non-Clinical Validation**
 - a. Analytical Measurement Range
 - i. Limits of Blank, Detection and Quantitation
 - ii. Linearity
 - b. Analytical Specificity
 - c. Precision
- V. **Clinical Validation**
 - a. Method Comparison with Predicate
 - b. Transference and Verification of Reference Interval (Venous)
 - c. Verification of Reference Interval with Finger Stick Samples
- VI. **Stability**
 - a. Reagent
 - b. Sample
 - c. Calibrators

Theranos Internal Only

Any retransmissions, reproductions, dissemination or other use of these materials by persons or entities other than the intended recipient is prohibited. This document supersedes all earlier or previous documents unless approved in writing.

TMP-00009 Rev. A, Released 08/01/13

theranos	LDT Validation Report	Theranos Phosphorus Assay	Rev:
		CL RPT-14057	1
Description	Validation Report for Modified Siemens Assay of Inorganic Phosphorus in Lithium Heparin Plasma		
Originator: Curtis Schneider		Date: 10/15/2013	

Overview

Eighty-eight percent of the phosphorus contained in the body is localized in bone in the form of hydroxyapatite. The remainder is involved in intermediary carbohydrate metabolism and in physiologically important substances such as phospholipids, nucleic acids, and adenosine triphosphate (ATP). Phosphorus occurs in blood in the form of inorganic phosphate and organically bound phosphoric acid. The small amount of extracellular organic phosphorus is found exclusively in the form of phospholipids. Serum contains approximately 2.5 to 4.5 mg/dL of inorganic phosphate (the fraction measure in routine biochemical assays). Serum phosphate concentrations are dependent on meals and variation in the secretion of hormones such as parathyroid hormone (PTH) and may vary widely.


Hypophosphatemia may have 4 general causes: shift of phosphate from extracellular to intracellular, renal phosphate wasting, loss from the gastrointestinal tract, and loss from intracellular stores.

Hyperphosphatemia is usually secondary to an inability of the kidneys to excrete phosphate. Other factors may relate to increased intake or a shift of phosphate from the tissues into the extracellular fluid.

I. Method Principle

Inorganic phosphorus reacts with ammonium molybdate in the presence of sulfuric acid to form an unreduced phosphomolybdate complex, which is measured as an endpoint reaction at 340/658 nm.

Reaction Equation



theranos	LDT Validation Report	Theranos Phosphorus Assay	Rev:
		CL RPT-14057	1
Description	Validation Report for Modified Siemens Assay of Inorganic Phosphorus in Lithium Heparin Plasma		
Originator: Curtis Schneider		Date: 10/15/2013	

II. Definitions and Abbreviations

The following definitions and abbreviations are used in this document and related documents and attachments:

- a. **Accuracy:** Accuracy is defined by CLSI as the closeness of agreement between a test result and an accepted reference value. Method accuracy is used in a different sense by the American Association of Pharmaceutical Scientists where it is expressed as percent relative error (%RE). Trueness, a related CLSI term, is the closeness of agreement between the average of a number of replicate measured quantity values and a reference quantity value.
- b. **Analyte:** Component represented in the name of a measurable quantity. The closely related term measure and is defined as the particular quantity subject to measurement.
- c. **Analytical sensitivity:** There are several alternative uses of this term. Most commonly, and for the purposes of this Validation Plan, it is used interchangeably with limit of detection. It is also used to describe the ability of an analytical method to assess small variations of the concentration of an analyte, such as the slope of the calibration curve (IUPAC).
- d. **Analytical specificity:** Ability of a test or procedure to correctly identify or quantify an entity, including in the presence of interfering substance(s) or phenomena.
- e. **Calibration:** Set of operations that establish, under specified conditions, the relationship between values of quantities indicated by a measuring instrument or measuring system, or values represented by a material measure or a reference material, and the corresponding values realized by standards. Under CLIA, calibration refers to the process of testing and adjusting an instrument, kit, or test system, to provide a known relationship between the measurement response and the value of the substance being measured by the test procedure (42 CFR 493.1217).
- f. **Calibrator:** A substance, material, or article intended to be used to establish the measurement relationships of a diagnostic medical device.
- g. **CLIA:** Clinical Laboratory Improvement Amendments of 1988. Congressional legislation that defined and requires specific quality assurance practices in clinical laboratories.
- h. **CLSI:** Clinical and Laboratory Standards Institute.

Theranos Confidential

Page 4 of 15

Any retransmissions, reproductions, dissemination or other use of these materials by persons or entities other than the intended recipient is prohibited. This document supersedes all earlier or previous documents unless approved in writing.

TMP-00009 Rev. A, Released 08/01/13

theranos	LDT Validation Report	Theranos Phosphorus Assay	Rev:
		CL RPT-14057	1
Description	Validation Report for Modified Siemens Assay of Inorganic Phosphorus in Lithium Heparin Plasma		
Originator: Curtis Schneider		Date: 10/15/2013	

- i. **Coefficient of Variation:** The ratio of the standard deviation to the average, often multiplied by 100 and expressed as a percentage, abbreviated as %CV .
- j. **Colorimetry:** A technique used to determine the concentration of colored compound(s) in solution.
- k. **Interfering substance:** A substance or quantity thereof that is not the measurand but that affects the result of the measurement.
- l. **IUPAC:** International Union of Pure and Applied Chemistry
- m. **LDT:** Laboratory –developed Test.
- n. **Linearity:** Linearity is the ability of a quantitative analytical method to provide results that are directly proportional to the concentrations of an analyte in test samples, within a given measuring interval. It is an important parameter to confirm when evaluating an analytical method because it verifies correct interpolation of results between points.
- o. **LMR:** Lower end of the measuring range is the lowest level at which defined conditions, including all stated characteristic of the method, are met.
- p. **LoB:** Limit of Blank is the highest value in a series of results on a sample that contains no analyte
- q. **LoD:** Limit of Detection is the lowest amount of analyte in a sample that can be detected with stated probability, although perhaps not quantified as an exact value.
- r. **LoQ:** When used without a prefix, the Limit of Quantitation is the lowest actual concentration at which an analyte is reliably detected and at which uncertainty of the test result is less than or equal to the goal set by the manufacturer or laboratory. The term may also be used with prefixes L for lower (LLOQ) and U for upper (ULOQ), respectively. Note: $LoB < LoD \leq LoQ$.
- s. **Matrix:** All components of a material system, except the analyte. A specimen matrix is the biological milieu in which an analyte exists (e.g., plasma, serum, urine, or other body fluids).

Any retransmissions, reproductions, dissemination or other use of these materials by persons or entities other than the intended recipient is prohibited. This document supersedes all earlier or previous documents unless approved in writing.

theranos	LDT Validation Report	Theranos Phosphorus Assay	Rev:
		CL RPT-14057	1
Description	Validation Report for Modified Siemens Assay of Inorganic Phosphorus in Lithium Heparin Plasma		
Originator: Curtis Schneider		Date: 10/15/2013	

- t. **Measuring Interval (reportable range; analytical measurement range or AMR):**
A measuring interval consists of all numeric values between the lower and upper numeric values for which a method can produce quantitative results suitable for clinical use. Where applicable, a linearity study is frequently used to establish or verify the measuring interval that can be reported for a measurement method. Alternatively, the lower limit of the measuring interval may be assigned as the LLOQ (LLOQ).
- u. **Precision:** Precision is the closeness of agreement between indications or measured quantity values obtained by replicate measurements on the same or similar objects under specified conditions. It is usually expressed numerically in terms of standard deviation (SD) or percent Coefficient of Variation (%CV).
- v. **Reference interval:** The interval between and including two reference limits. It is common practice to define a reference limit so a stated fraction of the reference values is less than or equal, or greater than or equal, to the respective upper or lower limit.
- w. **SOP:** Standard Operating Procedure.
- x. **Spectrophotometry:** The quantitative measurement of the transmission (or reflection) properties of a material as a function of wavelength.
- y. **Testing System:** The entirety of the testing process, including instrument, sample, reagents, supplies, and procedures. Personnel are sometimes included in the definition.

III. Pre-clinical Validation

a. Analytical Measurement Range

i. Limits of Blank, Detection and Quantitation

The limits of blank, detection, and quantitation were determined to be 0.0 mg/dL, 0.2 mg/dL and 1.0 mg/dL respectively.

theranos	LDT Validation Report	Theranos Phosphorus Assay	Rev:
		CL RPT-14057	1
Description	Validation Report for Modified Siemens Assay of Inorganic Phosphorus in Lithium Heparin Plasma		
Originator: Curtis Schneider		Date: 10/15/2013	

Limit of blank

CLSI guideline EP17-A section 4.3.1

Level	Number of samples	N	Mean	SD
Blank	1	20	0.00	0.00
Alpha Parametric LoB	5%		0.00	

Limit of detection

CLSI guideline EP17-A section 4.3.2

Level	Number of samples	N	Pooled SD
Low	1	20	0.12
Beta Parametric LoD	5%		0.20

Limit of quantitation

CLSI guideline EP17-A section 5.1

Level	Number of samples	N
Low	1	20
Bias	0.04	
Pooled imprecision	0.12	
95% total error	0.27	
Allowable error	-	

The lower limit of quantitation has been established at 1.0 mg/dL (11.4% CV and 105.4% recovery).

ii. Linearity

The Analytical Measurement Range (AMR) including linear measurement interval has been determined by Siemens. Refer to the Analytical Range section of the manufacturer product information insert for additional details.

Any retransmissions, reproductions, dissemination or other use of these materials by persons or entities other than the intended recipient is prohibited. This document supersedes all earlier or previous documents unless approved in writing.

TMP-00009 Rev. A, Released 08/01/13

theranos	LDT Validation Report	Theranos Phosphorus Assay	Rev:
		CL RPT-14057	1
Description	Validation Report for Modified Siemens Assay of Inorganic Phosphorus in Lithium Heparin Plasma		
Originator: Curtis Schneider		Date: 10/15/2013	

b. Analytical Specificity

The analytical specificity for this assay was determined by testing the effect of hemoglobin (100 mg/dL), bilirubin (10 mg/dL) and triglycerides (400 mg/dL) on plasma samples spiked with the interferents and then compared with un-spiked controls. Inorganic Phosphorus concentration at which the interference testing was performed at was 3 mg/dL. Non-interference was defined as the mean result from testing of spiked samples within 10% of the mean of the un-spiked samples. Recoveries were within 91.9% to 99.3% (see table below).

Table 1. Interference Testing For Inorganic Phosphorus

Analyte (mg/dL)	% Recovery		
	Interferent		
	Bilirubin (10 mg/dL)	Hemoglobin (100 mg/dL)	Triglycerides (400 mg/dL)
Phosphorus	98.0	99.5	91.9

No significant interference was observed.

c. Precision

theranos	LDT Validation Report	Theranos Phosphorus Assay	Rev:
		CL RPT-14057	1
Description	Validation Report for Modified Siemens Assay of Inorganic Phosphorus in Lithium Heparin Plasma		
Originator: Curtis Schneider		Date: 10/15/2013	

Level = Level 1

Number of observations	80
Number of runs	40
Number of days	20
Runs per day	2
Replicates per run	2
Mean	2.14

	SD	95% CI	CV	Allowable Total SD
Repeatability	0.07	0.06 to 0.09	3.2%	-
Between-run	0.00		0.0%	-
Between-day	0.03		1.2%	-
Within-laboratory	0.07	0.06 to 0.09	3.4%	0.21

Imprecision is less than allowable total imprecision: 10% upto 10mg/dL then 10%.

Level = Level 2

Number of observations	80
Number of runs	40
Number of days	20
Runs per day	2
Replicates per run	2
Mean	3.99

	SD	95% CI	CV	Allowable Total SD
Repeatability	0.10	0.09 to 0.13	2.6%	-
Between-run	0.00		0.0%	-
Between-day	0.09		2.3%	-
Within-laboratory	0.14	0.12 to 0.17	3.5%	0.40

Imprecision is less than allowable total imprecision: 10% upto 10mg/dL then 10%.

Theranos Confidential

Page 9 of 15

Any retransmissions, reproductions, dissemination or other use of these materials by persons or entities other than the intended recipient is prohibited. This document supersedes all earlier or previous documents unless approved in writing.

TMP-00009 Rev. A, Released 08/01/13

theranos	LDT Validation Report	Theranos Phosphorus Assay	Rev:
		CL RPT-14057	1
Description	Validation Report for Modified Siemens Assay of Inorganic Phosphorus in Lithium Heparin Plasma		
Originator: Curtis Schneider		Date: 10/15/2013	

Level = Level 3

Number of observations	80
Number of runs	40
Number of days	20
Runs per day	2
Replicates per run	2
Mean	7.28

	SD	95% CI	CV	Allowable Total SD
Repeatability	0.14	0.12 to 0.18	2.0%	-
Between-run	0.04		0.6%	-
Between-day	0.12		1.7%	-
Within-laboratory	0.19	0.16 to 0.24	2.6%	0.73

Imprecision is less than allowable total imprecision: 10% upto 10mg/dL then 10%.

The percent CV reported as zeros in the above precision summary are most likely a consequence of rounding the values in StatisPro.

IV. Clinical Validation

a. Method Comparison with Predicate (Accuracy/Comparability)

To test the accuracy of the assay on the Theranos System, forty nine (49) unique patient samples were screened on the predicate method (Siemens, Advia) and on the Theranos method. Using the predicate method twenty eight (28) values were within the reference range (2.4 — 5.1 mg/dL), one (1) was below the reference range, and twenty (20) were above the reference range. Based on the results of the data examination, either a simple linear regression or alternative procedures were used to estimate expected (average) bias and the confidence interval of expected bias at the desired medical decision level(s) as per CLSI guidance EP09-A2. StatisPro was used for bias calculations. These estimates were compared with internal criteria to judge the acceptability of the Theranos method. Each sample was run in duplicate on the predicate, and the average used for comparison to the Theranos method. Some samples were stored before analysis on both methods. If the confidence interval for the predicted bias includes the defined acceptable bias or if the acceptable bias is

Theranos Confidential

Page 10 of 15

Any retransmissions, reproductions, dissemination or other use of these materials by persons or entities other than the intended recipient is prohibited. This document supersedes all earlier or previous documents unless approved in writing.

TMP-00009 Rev. A, Released 08/01/13

theranos	LDT Validation Report	Theranos Phosphorus Assay	Rev:
		CL RPT-14057	1
Description	Validation Report for Modified Siemens Assay of Inorganic Phosphorus in Lithium Heparin Plasma		
Originator: Curtis Schneider		Date: 10/15/2013	

greater than the higher limit of the confidence interval of the predicted bias, then the data do not show that the bias of the Theranos method is different from the acceptable bias or there is a high probability (97%) that the predicated bias is acceptable, respectively. The acceptable bias at each medical decision level was determined based on the total allowable error (TEa) minus the measured precision at the level closest to that decision level. Total allowable error (TEa) was taken from American Proficiency Institute (API) peer proficiency testing criteria or CLIA proficiency testing criteria for acceptable analytical performance, as printed in the Federal Register February 28, 1992;57(40):7002-186, when available. The TEa for Phosphorus is 10%. The table below shows the allowable bias and precision at 3 levels (values shown in parentheses) and the corresponding closest medical decision limits.

Table 2. Allowable Bias and Precision at the Medical Decision Levels

Medical Decision Levels (mg/dL)	1.5(2.1)	2.5(4.0)	5.0(7.3)
Precision (%)	1.2	2.3	1.7
Allowable Bias (%)	8.8	7.7	8.3

theranos	LDT Validation Report	Theranos Phosphorus Assay	Rev:
		CL RPT-14057	1
Description	Validation Report for Modified Siemens Assay of Inorganic Phosphorus in Lithium Heparin Plasma		
Originator: Curtis Schneider		Date: 10/15/2013	




Figure 1. Graph showing Theranos method versus Predicate Method (Siemens Advia).

Simple linear regression was used to establish a slope, intercept and an r^2 . The slope, intercept and clinical correlation were determined to be 0.97, 0.12 and 0.99 respectively.

Comparability

(CLSI guideline EP09-A3-IR section 7)

Level ID	Value	Difference	SE	95% CI	Allowable difference
1.500000	0.0700872	0.05525847	-0.0416685 to 0.1812629	0.1155000	
2.500000	0.0390539	0.04319135	-0.0478359 to 0.1256437	0.2250000	
5.000000	-0.0385543	0.03059695	-0.1001074 to 0.0229986	0.4500000	

Difference is less than allowable bias: 7.7% upto 2.5mg/dL then 9%.

theranos	LDT Validation Report	Theranos Phosphorus Assay	Rev:
		CL RPT-14057	1
Description	Validation Report for Modified Siemens Assay of Inorganic Phosphorus in Lithium Heparin Plasma		
Originator: Curtis Schneider	Date: 10/15/2013		

The difference between the two methods is not greater than the allowable difference. The performance requirement is verified.

b. Transference and Verification of Reference Interval (Venous)

Reference ranges were modified by applying the regression equation to the lower and upper reference limits of existing reference interval to generate a new reference range. New reference ranges were verified using a minimum of twenty (20) new normal subjects

New reference ranges were verified using a total of fifty (50) new normal subjects with matched Lithium heparin venous and finger sticks samples. For a reference range to pass verification, 95% of values should fall within the upper and lower reference limits and 5% or fewer values fall outside of the upper and lower reference limits. For venous verification 50 (100%) values fell within the new reference range and 0 (0%) values fell outside the new reference range. See graphs below for venous samples verification.




Figure 2. Graph showing venous sample reference range verification.

theranos	LDT Validation Report	Theranos Phosphorus Assay	Rev:
		CL RPT-14057	1
Description	Validation Report for Modified Siemens Assay of Inorganic Phosphorus in Lithium Heparin Plasma		
Originator: Curtis Schneider		Date: 10/15/2013	

c. Verification of Reference Interval with Finger Stick Samples

New reference ranges were also verified with venous matched finger sticks (Lithium heparin) samples from a total of forty six (46) new normal subjects. The finger stick samples were collected in a Theranos blood collection device (BCD) configured with separate Lithium heparin and EDTA vessels. For finger stick verification 45 values (98%) fell within the new reference range and 1 value (2%) fell outside the new reference range. See graphs below for finger stick samples verification.




Figure 3. Graph showing Finger stick sample reference range verification.

The new reference range for finger stick inorganic phosphorus was determined to be 2.5 – 5.1 mg/dL.

VI. Stability

a. Reagents

On-board Reagent Stability

Any retransmissions, reproductions, dissemination or other use of these materials by persons or entities other than the intended recipient is prohibited. This document supersedes all earlier or previous documents unless approved in writing.

theranos	LDT Validation Report	Theranos Phosphorus Assay	Rev:
		CL RPT-14057	1
Description	Validation Report for Modified Siemens Assay of Inorganic Phosphorus in Lithium Heparin Plasma		
Originator: Curtis Schneider	Date: 10/15/2013		

System	Stability
ADVIA 1200	50 days
ADVIA 1650/1800	50 days
ADVIA 2400	60 days

For all systems, unopened reagents are stable until the expiration date printed on the product label when stored at 15°C - 25°C. Do not freeze the reagents.

For additional details, refer to the Methods Introduction section of the system-specific Operator's Guide.

b. Sample

Plasma samples for inorganic phosphorous analysis are stable for 2 weeks at 2-8 °C, or at least 90 days at -20 °C.

c. Calibrators

Siemens Chemistry Calibrators should be stored at 2-8 °C, protected from light, and are stable until the expiration date on the vial label. Opened calibrators are stable for 48 hours, except for total and direct bilirubin, which are stable for 8 hours.

REVISION HISTORY			
Revision Level	Effective Date	Initiator	ECO Number
A	11/10/2013	A. Rosendorff	CL ECO-00118
Section Number	Description and Justification of Changes		
All	Initial Release		

Any retransmissions, reproductions, dissemination or other use of these materials by persons or entities other than the intended recipient is prohibited. This document supersedes all earlier or previous documents unless approved in writing.

TMP-00009 Rev. A, Released 08/01/13