


<b>theranos</b>	LDT Validation Report	Theranos Lactate Assay	Rev:
		<b>CL-RPT-14071</b>	<b>1</b>
Description	Validation Report for Modified Siemens Assay of Lactate in Lithium Heparin Plasma		
Originator: Curtis Schneider		Date: 10/15/2013	


### Validation of Modified Siemens Lactate Assay

**Author(s):**


Signature: 	Date: 12/12/13
Name: Paul Patel, Ph.D.	Title: Team Lead, General Chemistry

**Reviewer(s):**

Signature:	Date:
Name:	Title:

Signature: 	Date: 1/25/2014
Name: Daniel Young, Ph.D.	Title: Vice President

**Approver(s):**

Signature: 	Date: 12/12/2013
Name: Adam Rosendorff, M.D.	Title: Laboratory Director

 9/19/13

Sunil S. Dhawan M.D.

**Theranos Confidential**
**Page 1 of 16**

Any retransmissions, reproductions, dissemination or other use of these materials by persons or entities other than the intended recipient is prohibited. This document supersedes all earlier or previous documents unless approved in writing.

TMP-00009 Rev. A, Released 08/01/13

<b>theranos</b>	<b>LDT Validation Report</b>	Theranos Lactate Assay	Rev:
		<b>CL-RPT-14071</b>	<b>1</b>
Description	Validation Report for Modified Siemens Assay of Lactate in Lithium Heparin Plasma		
Originator: Curtis Schneider		Date: 10/15/2013	

## Lactate Plasma Assay

- I. Overview**
- II. Method Principle**
- III. Definitions and Abbreviations**
- IV. Non-Clinical Validation**
  - a. Analytical Measurement Range
    - i. Limits of Blank, Detection and Quantitation
    - ii. Linearity
  - b. Analytical Specificity
  - c. Precision
- V. Clinical Validation**
  - a. Method Comparison with Predicate
  - b. Transference and Verification of Reference Interval (Venous)
  - c. Verification of Reference Interval with Finger Stick Samples
- VI. Stability**
  - a. Reagent
  - b. Sample
  - c. Calibrators

**Theranos Confidential**

**Page 2 of 16**

Any retransmissions, reproductions, dissemination or other use of these materials by persons or entities other than the intended recipient is prohibited. This document supersedes all earlier or previous documents unless approved in writing.

TMP-00009 Rev. A, Released 08/01/13

<b>theranos</b>	<b>LDT Validation Report</b>	Theranos Lactate Assay	Rev:
		<b>CL-RPT-14071</b>	<b>1</b>
Description	Validation Report for Modified Siemens Assay of Lactate in Lithium Heparin Plasma		
Originator: Curtis Schneider		Date: 10/15/2013	

## Overview


Lactate is the end product of anaerobic carbohydrate metabolism. Major sites of production are skeletal muscle, brain, and erythrocytes. Lactate is metabolized by the liver. The concentration of lactate depends on the rate of production and the rate of liver clearance. The liver can adequately clear lactate until the concentration reaches approximately 2 mmol/L. When this level is exceeded, lactate begins to accumulate rapidly. For example, while resting lactate levels are usually <1 mmol/L, during strenuous exercise, levels can rise >20 mmol/L within a few seconds.

Lactic acidosis signals the deterioration of the cellular oxidative process and is associated with hyperpnea, weakness, fatigue, stupor, and finally coma. These conditions may be irreversible, even after treatment is administered. Lactate acidosis may be associated with hypoxic conditions (eg, shock, hypovolemia, heart failure, pulmonary insufficiency), metabolic disorders (eg, diabetic ketoacidosis, malignancies), and toxin exposures (eg, ethanol, methanol, salicylates).

## I. Method Principle

Lactate is oxidized by lactate oxidase to pyruvate and hydrogen peroxide. Lactate is measured by the formation of dye from hydrogen peroxide and a chromogen in presence of a peroxidase. The corresponding change in absorbance at 545/694 nm is proportional to the plasma lactate concentration.

### Reaction Equation



\*TOOS = N-Ethyl-N-(2-Hydroxy-3-sulphopropyl) M-Toluidine

<b>theranos</b>	LDT Validation Report	Theranos Lactate Assay	Rev:
		<b>CL-RPT-14071</b>	<b>1</b>
Description	Validation Report for Modified Siemens Assay of Lactate in Lithium Heparin Plasma		
Originator: Curtis Schneider		Date: 10/15/2013	

## II. Definitions and Abbreviations

The following definitions and abbreviations are used in this document and related documents and attachments:

- a. **Accuracy:** Accuracy is defined by CLSI as the closeness of agreement between a test result and an accepted reference value. Method accuracy is used in a different sense by the American Association of Pharmaceutical Scientists where it is expressed as percent relative error (%RE). Trueness, a related CLSI term, is the closeness of agreement between the average of a number of replicate measured quantity values and a reference quantity value.
- b. **Analyte:** Component represented in the name of a measurable quantity. The closely related term measureand is defined as the particular quantity subject to measurement.
- c. **Analytical sensitivity:** There are several alternative uses of this term. Most commonly, and for the purposes of this Validation Plan, it is used interchangeably with limit of detection. It is also used to describe the ability of an analytical method to assess small variations of the concentration of an analyte, such as the slope of the calibration curve (IUPAC).
- d. **Analytical specificity:** Ability of a test or procedure to correctly identify or quantify an entity, including in the presence of interfering substance(s) or phenomena.
- e. **Calibration:** Set of operations that establish, under specified conditions, the relationship between values of quantities indicated by a measuring instrument or measuring system, or values represented by a material measure or a reference material, and the corresponding values realized by standards. Under CLIA, calibration refers to the process of testing and adjusting an instrument, kit, or test system, to provide a known relationship between the measurement response and the value of the substance being measured by the test procedure (42 CFR 493.1217).
- f. **Calibrator:** A substance, material, or article intended to be used to establish the measurement relationships of a diagnostic medical device.
- g. **CLIA:** Clinical Laboratory Improvement Amendments of 1988. Congressional legislation that defined and requires specific quality assurance practices in clinical laboratories.
- h. **CLSI:** Clinical and Laboratory Standards Institute.

**Theranos Confidential**

**Page 4 of 16**

Any retransmissions, reproductions, dissemination or other use of these materials by persons or entities other than the intended recipient is prohibited. This document supersedes all earlier or previous documents unless approved in writing.

TMP-00009 Rev. A, Released 08/01/13

<b>theranos</b>	LDT Validation Report	Theranos Lactate Assay	Rev:
		<b>CL-RPT-14071</b>	<b>1</b>
Description	Validation Report for Modified Siemens Assay of Lactate in Lithium Heparin Plasma		
Originator: Curtis Schneider		Date: 10/15/2013	

- i. **Coefficient of Variation:** The ratio of the standard deviation to the average, often multiplied by 100 and expressed as a percentage, abbreviated as %CV .
- j. **Colorimetry:** A technique used to determine the concentration of colored compound(s) in solution.
- k. **Interfering substance:** A substance or quantity thereof that is not the measurand but that affects the result of the measurement.
- l. **IUPAC:** International Union of Pure and Applied Chemistry
- m. **LDT:** Laboratory –developed Test.
- n. **Linearity:** Linearity is the ability of a quantitative analytical method to provide results that are directly proportional to the concentrations of an analyte in test samples, within a given measuring interval. It is an important parameter to confirm when evaluating an analytical method because it verifies correct interpolation of results between points.
- o. **LMR:** Lower end of the measuring range is the lowest level at which defined conditions, including all stated characteristic of the method, are met.
- p. **LoB:** Limit of Blank is the highest value in a series of results on a sample that contains no analyte.
- q. **LoD:** Limit of Detection is the lowest amount of analyte in a sample that can be detected with stated probability, although perhaps not quantified as an exact value.
- r. **LoQ:** When used without a prefix, the Limit of Quantitation is the lowest actual concentration at which an analyte is reliably detected and at which uncertainty of the test result is less than or equal to the goal set by the manufacturer or laboratory. The term may also be used with prefixes L for lower (LLOQ) and U for upper (ULOQ), respectively. Note:  $LoB < LoD \leq LoQ$ .
- s. **Matrix:** All components of a material system, except the analyte. A specimen matrix is the biological milieu in which an analyte exists (e.g., plasma, serum, urine, or other body fluids).

**Theranos Confidential**

**Page 5 of 16**

Any retransmissions, reproductions, dissemination or other use of these materials by persons or entities other than the intended recipient is prohibited. This document supersedes all earlier or previous documents unless approved in writing.

TMP-00009 Rev. A, Released 08/01/13

<b>theranos</b>	<b>LDT Validation Report</b>	Theranos Lactate Assay	Rev:
		<b>CL-RPT-14071</b>	<b>1</b>
Description	Validation Report for Modified Siemens Assay of Lactate in Lithium Heparin Plasma		
Originator: Curtis Schneider		Date: 10/15/2013	

- t. **Measuring Interval (reportable range; analytical measurement range or AMR):**  
A measuring interval consists of all numeric values between the lower and upper numeric values for which a method can produce quantitative results suitable for clinical use. Where applicable, a linearity study is frequently used to establish or verify the measuring interval that can be reported for a measurement method. Alternatively, the lower limit of the measuring interval may be assigned as the LLOQ (LLOQ).
- u. **Precision:** Precision is the closeness of agreement between indications or measured quantity values obtained by replicate measurements on the same or similar objects under specified conditions. It is usually expressed numerically in terms of standard deviation (SD) or percent Coefficient of Variation (%CV).
- v. **Reference interval:** The interval between and including two reference limits. It is common practice to define a reference limit so a stated fraction of the reference values is less than or equal, or greater than or equal, to the respective upper or lower limit.
- w. **SOP:** Standard Operating Procedure.
- x. **Spectrophotometry:** The quantitative measurement of the transmission (or reflection) properties of a material as a function of wavelength.
- y. **Testing System:** The entirety of the testing process, including instrument, sample, reagents, supplies, and procedures. Personnel are sometimes included in the definition.

### III. Pre-clinical Validation

#### a. Analytical Measurement Range

##### i. Limits of Blank, Detection and Quantitation

The limits of blank, detection, and quantitation were determined to be 0.1 mmol/L, 0.13 mmol/L and 0.13 mmol/L respectively.

<b>theranos</b>	LDT Validation Report	Theranos Lactate Assay	Rev:
		<b>CL-RPT-14071</b>	<b>1</b>
Description	Validation Report for Modified Siemens Assay of Lactate in Lithium Heparin Plasma		
Originator: Curtis Schneider		Date: 10/15/2013	

### Limit of blank

CLSI guideline EP17-A section 4.3.1

Level	Number of samples	N	Mean	SD
Blank	1	20	0.02	0.05
Alpha	5%			
Parametric LoB	0.10			

### Limit of detection

CLSI guideline EP17-A section 4.3.2

Level	Number of samples	N	Pooled SD
Low	1	20	0.02
Beta	5%		
Parametric LoD	0.13		

### Limit of quantitation

CLSI guideline EP17-A section 5.1

Level	Number of samples	N
Low	1	20
Bias	-0.11	
Pooled imprecision	0.02	
95% total error	-0.15	
Allowable error	0.4	
LoQ	0.13	

95% total error is less than allowable error: 20%.

LoQ has been established.

**Theranos Confidential**

**Page 7 of 16**


Any retransmissions, reproductions, dissemination or other use of these materials by persons or entities other than the intended recipient is prohibited. This document supersedes all earlier or previous documents unless approved in writing.

TMP-00009 Rev. A, Released 08/01/13

<b>theranos</b>	LDT Validation Report	Theranos Lactate Assay	Rev:
		<b>CL-RPT-14071</b>	<b>1</b>
Description	Validation Report for Modified Siemens Assay of Lactate in Lithium Heparin Plasma		
Originator: Curtis Schneider		Date: 10/15/2013	

## ii. Linearity

The Analytical Measurement Range (AMR) including linear measurement interval has been determined for Lactate in plasma. This method is linear from 0.2 – 10.9 mM (1.6 – 98.6 mg/dL) within the 10% allowable non-linearity in this interval.



A linear relationship fits the data better than a nonlinear relationship over the measuring interval.

Level	Mean	Linear fit	Nonlinear fit	Nonlinearity	Allowable nonlinearity
1	1.60	1.60	-	-	0.16
2	23.90	23.90	-	-	2.39
3	48.30	48.30	-	-	4.83
4	72.90	72.90	-	-	7.29
5	98.60	98.60	-	-	9.86

A linear relationship fits the data better than a nonlinear relationship over the measuring interval.

Nonlinearity is less than allowable nonlinearity: 10%.  
Performance requirement verified over the measuring interval.

**Theranos Confidential**

**Page 8 of 16**

Any retransmissions, reproductions, dissemination or other use of these materials by persons or entities other than the intended recipient is prohibited. This document supersedes all earlier or previous documents unless approved in writing.

TMP-00009 Rev. A, Released 08/01/13



<b>theranos</b>	<b>LDT Validation Report</b>	Theranos Lactate Assay	Rev:
		<b>CL-RPT-14071</b>	<b>1</b>
Description	Validation Report for Modified Siemens Assay of Lactate in Lithium Heparin Plasma		
Originator: Curtis Schneider		Date: 10/15/2013	

### Analytical Specificity

The analytical specificity for this assay was determined by testing the effect of hemoglobin (100 mg/dL), bilirubin (10 mg/dL) and triglycerides (400 mg/dL) on plasma samples spiked with the interferents and then compared with un-spiked controls. Lactate concentration at which the interference testing was performed at was 24.7 mmol/L (222 mg/dL). Non-interference was defined as the mean result from testing of spiked samples within 10% of the mean of the un-spiked samples. Recoveries were within 97.4% to 99.3% (see table below).

Table 1. Interference Testing For Lactate.

Analyte (mg/dL)	% Recovery		
	Interferent		
	Bilirubin (10 mg/dL)	Hemoglobin (100 mg/dL)	Triglycerides (400 mg/dL)
Lactate	99.3	99.0	97.4

No significant interference was observed.

### b. Precision

Level = Level 1

Number of observations	80
Number of runs	40
Number of days	20
Runs per day	2
Replicates per run	2
Mean	8.55

	SD	95% CI	CV	Allowable Total SD
Repeatability	0.42	0.35 to 0.54	5.0%	-
Between-run	0.13		1.5%	-
Between-day	0.52		6.0%	-
Within-laboratory	0.68	0.56 to 0.88	8.0%	0.85

Imprecision is less than allowable total imprecision: 10% upto 20mg/dL then 10%.

**Theranos Confidential**

**Page 9 of 16**

Any retransmissions, reproductions, dissemination or other use of these materials by persons or entities other than the intended recipient is prohibited. This document supersedes all earlier or previous documents unless approved in writing.

TMP-00009 Rev. A, Released 08/01/13

<b>theranos</b>	<b>LDT Validation Report</b>	Theranos Lactate Assay	Rev:
		<b>CL-RPT-14071</b>	<b>1</b>
Description	Validation Report for Modified Siemens Assay of Lactate in Lithium Heparin Plasma		
Originator: Curtis Schneider		Date: 10/15/2013	

**Level = Level 2**

Number of observations 80  
Number of runs 40  
Number of days 20  
Runs per day 2  
Replicates per run 2

Mean 25.86

	SD	95% CI	CV	Allowable Total SD
Repeatability	0.74	0.61 to 0.95	2.9%	-
Between-run	0.00		0.0%	-
Between-day	1.84		7.1%	-
Within-laboratory	1.99	1.54 to 2.79	7.7%	2.59

Imprecision is less than allowable total imprecision: 10% upto 20mg/dL then 10%.

**Level = Level 3**

Number of observations 76  
Number of runs 38  
Number of runs excluded 2  
Number of days 20  
% of days with 1 run 10%  
Runs per day 2  
Replicates per run 2

**Level = Level 3**

Number of observations 76  
Number of runs 38  
Number of runs excluded 2  
Number of days 20  
% of days with 1 run 10%  
Runs per day 2  
Replicates per run 2

CLSI guideline EP05-A2 section 10.4 recommends a minimum of 40 runs, with 2 replicates per run.

Mean 47.24

	SD	95% CI	CV	Allowable Total SD
Repeatability	0.98	0.81 to 1.26	2.1%	-
Between-run	0.22		0.5%	-
Between-day	3.29		7.0%	-
Within-laboratory	3.44	2.65 to 4.88	7.3%	4.72

Imprecision is less than allowable total imprecision: 10% upto 20mg/dL then 10%.

**Theranos Confidential**

**Page 10 of 16**

Any retransmissions, reproductions, dissemination or other use of these materials by persons or entities other than the intended recipient is prohibited. This document supersedes all earlier or previous documents unless approved in writing.

TMP-00009 Rev. A, Released 08/01/13

<b>theranos</b>	LDT Validation Report	Theranos Lactate Assay	Rev:
		<b>CL-RPT-14071</b>	<b>1</b>
Description	Validation Report for Modified Siemens Assay of Lactate in Lithium Heparin Plasma		
Originator: Curtis Schneider		Date: 10/15/2013	

The percent CV reported as zeros in the above precision summary are most likely a consequence of rounding the values in StatisPro.

#### IV. Clinical Validation

##### a. Method Comparison with Predicate (Accuracy/Comparability)

To test the accuracy of the assay on the Theranos System, forty-six (46) unique patient samples were screened on the predicate method (Siemens, Advia) and on the Theranos method. Eight samples were excluded as outliers (mean absolute difference greater than 4). Using the predicate method twenty (20) values were within the reference range (0.5 – 2.2 mM (4.5 -19.8 mg/dL)) and eighteen (18) were above the reference range. Based on the results of the data examination, either a simple linear regression or alternative procedures were used to estimate expected (average) bias and the confidence interval of expected bias at the desired medical decision level(s) as per CLSI guidance EP09-A2. StatisPro was used for bias calculations. These estimates were compared with internal criteria to judge the acceptability of the Theranos method. Each sample was run in duplicate on the predicate, and the average used for comparison to the Theranos method. Some samples were stored before analysis on both methods. If the confidence interval for the predicted bias includes the defined acceptable bias or if the acceptable bias is greater than the higher limit of the confidence interval of the predicted bias, then the data do not show that the bias of the Theranos method is different from the acceptable bias or there is a high probability (97%) that the predicated bias is acceptable, respectively. The acceptable bias at each medical decision level was determined based on the total allowable error (TEa) minus the measured precision at the level closest to that decision level. Total allowable error (TEa) was taken from American Proficiency Institute (API) peer proficiency testing criteria or CLIA proficiency testing criteria for acceptable analytical performance, as printed in the Federal Register February 28, 1992;57(40):7002-186, when available. The TEa for Lactate is 10%. The table below shows the allowable bias and precision at 2 level (actual values shown in parentheses) and the corresponding closest medical decision limits.

**Theranos Confidential**

**Page 11 of 16**


Any retransmissions, reproductions, dissemination or other use of these materials by persons or entities other than the intended recipient is prohibited. This document supersedes all earlier or previous documents unless approved in writing.

TMP-00009 Rev. A, Released 08/01/13

<b>theranos</b>	LDT Validation Report	Theranos Lactate Assay	Rev:
		<b>CL-RPT-14071</b>	<b>1</b>
Description	Validation Report for Modified Siemens Assay of Lactate in Lithium Heparin Plasma		
Originator: Curtis Schneider		Date: 10/15/2013	

Table 2. Allowable Bias and Precision at the Medical Decision Levels

Medical Decision Levels ( mmol/L)	2.2 (2.9)	5.6 (5.2)
Precision (%)	7.1	7.0
Allowable Bias (%)	2.9	3.0



**Figure 1.** Graph showing Theranos method versus Predicate Method (Siemens Advia).

**Theranos Confidential**

**Page 12 of 16**

Any retransmissions, reproductions, dissemination or other use of these materials by persons or entities other than the intended recipient is prohibited. This document supersedes all earlier or previous documents unless approved in writing.

TMP-00009 Rev. A, Released 08/01/13

<b>theranos</b>	LDT Validation Report	Theranos Lactate Assay	Rev:
		<b>CL-RPT-14071</b>	<b>1</b>
Description	Validation Report for Modified Siemens Assay of Lactate in Lithium Heparin Plasma		
Originator: Curtis Schneider		Date: 10/15/2013	

Simple linear regression was used to establish a slope, intercept and an  $r^2$ . The slope, intercept and clinical correlation were determined to be 1.07, -0.24 and 0.99 respectively.

#### Comparability

CLSI guideline EP06-A2-IR section 7

Level ID	Value	Difference	SE	95% CI	Allowable difference
	5.000000	0.1263119	0.06129991	0.0019900 to 0.2506339	0.5000000


Difference is less than allowable bias: 10% upto 10mmol/L then 10%.

The difference between the two methods is not greater than the allowable difference. The performance requirement is verified.

#### b. Establishment of Reference Interval with Venous Samples

A new reference range was established for Lactate in Lithium Heparin plasma according to CLSI guideline C28-A3 by testing 120 normal subjects (67 males, 71 females).

Frequency histogram  
CLSI guideline C28-A3 section 8.2




**Figure 2.** Frequency histogram showing distribution of venous values for reference range establishment.

<b>theranos</b>	<b>LDT Validation Report</b>	Theranos Lactate Assay	Rev:
		<b>CL-RPT-14071</b>	<b>1</b>
<b>Description</b>	Validation Report for Modified Siemens Assay of Lactate in Lithium Heparin Plasma		
<b>Originator: Curtis Schneider</b>		<b>Date: 10/15/2013</b>	

#### Reference Interval

CLSI guideline C28-A3 section 9.4.1/9.5.1



Computation method	95% reference interval	Lower reference limit 90% CI	Upper reference limit 90% CI
Nonparametric	5.00 to 20.47	4.20 to 5.30	18.60 to 22.10


The new expected venous reference values for this method are 0.55 - 2.28 mM (5.0 – 20.5 mg/dL).

#### c. Establishment of Reference Interval with Finger Stick Samples

A new reference range was also established for Lactate in Lithium Heparin plasma according to CLSI guideline C28-A3 by testing matched finger samples from the same 120 normal subjects (67 males, 71 females) mentioned above.

#### Frequency histogram

CLSI guideline C28-A3 section 9.2




**Figure 3.** Frequency histogram showing distribution of finger stick values for reference range establishment.

<b>theranos</b>	LDT Validation Report	Theranos Lactate Assay	Rev:
		<b>CL-RPT-14071</b>	<b>1</b>
Description	Validation Report for Modified Siemens Assay of Lactate in Lithium Heparin Plasma		
Originator: Curtis Schneider		Date: 10/15/2013	

#### Reference Interval

CLSI guideline C28-A3 section 9.4.1/9.5.1



Computation method	95% reference interval	Lower reference limit 90% CI	Upper reference limit 90% CI
Nonparametric	6.19 to 31.58	5.40 to 7.70	27.20 to 33.00

The new expected reference range values for finger stick lactate was determined to be 0.7 – 3.5 mM (6.2 – 31.6 mg/dL).

## VI. Stability

### a. Reagents

#### On-board Reagent Stability

System	Stability
ADVIA 1200	60 days
ADVIA 1650/1800	45 days
ADVIA 2400	45 days

**Theranos Confidential**

**Page 15 of 16**

Any retransmissions, reproductions, dissemination or other use of these materials by persons or entities other than the intended recipient is prohibited. This document supersedes all earlier or previous documents unless approved in writing.

TMP-00009 Rev. A, Released 08/01/13

<b>theranos</b>	LDT Validation Report	Theranos Lactate Assay	Rev:
		<b>CL-RPT-14071</b>	<b>1</b>
Description	Validation Report for Modified Siemens Assay of Lactate in Lithium Heparin Plasma		
Originator: Curtis Schneider		Date: 10/15/2013	

For all systems, unopened reagents are stable until the expiration date printed on the product label when stored at 2°C - 8°C and protected from light with screw caps tightly closed. Do not freeze the reagents.

For additional details, refer to the Methods Introduction section of the system-specific Operator's Guide.

**b. Sample**

Plasma samples for lactate analysis are stable for 1 week at 2-8 °C, or at least 2 week at -20 °C.

**c. Calibrators**

Siemens Special Chemistry Calibrators should be stored at 2-8 °C, protected from light, and are stable until the expiration date on the vial label. Opened reconstituted calibrators are stable for 7 days, except for acid phosphatase, which is stable for at least 2 days.

Theranos Internal Only