


<b>theranos</b>	LDT Validation Report	Theranos Total Protein Assay	Rev:
		<b>CL-RPT-14038</b>	<b>1</b>
Description	Validation Report for Modified Siemens Assay of Total Protein in Lithium Heparin Plasma		
Originator: Curtis Schneider		Date: 09/24/2013	


### Validation of Modified Siemens Total Protein Assay

#### Author(s):


Signature: 	Date: 11/17/13
Name: Paul Patel, Ph.D.	Title: Team Lead, General Chemistry

#### Reviewer(s):

Signature:	Date:
Name:	Title:

Signature:  SURAJ SAKSENA	Date:
Name: Daniel Young, Ph.D.	Title: Vice President

#### Approver(s):

Signature: 	Date: 11/7/13
Name: Adam Rosendorff, M.D.	Title: Laboratory Director

 9/19/13

Sunil S. Dhawan M.D.

**Theranos Confidential**

Page 1 of 15

Any retransmissions, reproductions, dissemination or other use of these materials by persons or entities other than the intended recipient is prohibited. This document supersedes all earlier or previous documents unless approved in writing.

TMP-00009 Rev. A, Released 08/01/13

<b>theranos</b>	LDT Validation Report	Theranos Total Protein Assay <b>CL-RPT-14038</b>	Rev: <b>1</b>
	Description	Validation Report for Modified Siemens Assay of Total Protein in Lithium Heparin Plasma	
Originator: Curtis Schneider		Date: 09/24/2013	

## Total Protein Plasma Assay

- I. Overview**
- II. Method Principle**
- III. Definitions and Abbreviations**
- IV. Pre-clinical Validation**
  - a. Analytical Measurement Range
    - i. Limits of Blank, Detection and Quantitation
    - ii. Linearity
  - b. Analytical Specificity
  - c. Precision
- V. Clinical Validation**
  - a. Method Comparison with Predicate
  - b. Transference and Verification of Reference Interval (Venous)
  - c. Verification of Reference Interval with Finger Stick Samples
- VI. Stability**
  - a. Reagent
  - b. Sample
  - c. Calibrators

Theranos Internal Only

<b>theranos</b>	LDT Validation Report	Theranos Total Protein Assay <b>CL-RPT-14038</b>	Rev: <b>1</b>
	Validation Report for Modified Siemens Assay of Total Protein in Lithium Heparin Plasma		
Description		Validation Report for Modified Siemens Assay of Total Protein in Lithium Heparin Plasma	
Originator: Curtis Schneider		Date: 09/24/2013	

## I. Overview

Plasma proteins are synthesized predominantly in the liver; immunoglobulins are synthesized by mononuclear cells of lymph nodes, spleen and bone marrow. The two general causes for alterations of serum total protein concentrations are a change in the volume of plasma water and/or a change in the concentration of one or more of the specific proteins in the plasma. Of the individual serum proteins, albumin is present in such high concentrations that low levels of this protein alone may cause hypoproteinemia.

Hemoconcentration (decrease in the volume of plasma water) results in relative hyperproteinemia whereas hemodilution results in relative hypoproteinemia. In both situations, the relative concentrations of all the individual plasma proteins are affected to the same degree.

Hyperproteinemia may be seen during dehydration, due to inadequate water intake, to excessive water loss (e.g., severe vomiting, diarrhea, Addison's disease and diabetic acidosis), or as a result of increased production of proteins. Increased polyclonal protein production is seen in reactive, inflammatory processes and increased monoclonal protein production is seen in some hematopoietic neoplasms (e.g., multiple myeloma, Waldenstrom's macroglobulinemia, monoclonal gammopathy of undetermined significance).

## II. Method Principle

Protein peptide bonds interact with the cupric ions to form a purple complex that is measured as an endpoint reaction at 545 nm.

### Reaction Equation



<b>theranos</b>	LDT Validation Report	Theranos Total Protein Assay	Rev:
		<b>CL-RPT-14038</b>	<b>1</b>
Description	Validation Report for Modified Siemens Assay of Total Protein in Lithium Heparin Plasma		
Originator: Curtis Schneider		Date: 09/24/2013	

### III. Definitions and Abbreviations

The following definitions and abbreviations are used in this document and related documents and attachments:

- a. **Accuracy:** Accuracy is defined by CLSI as the closeness of agreement between a test result and an accepted reference value. Method accuracy is used in a different sense by the American Association of Pharmaceutical Scientists where it is expressed as percent relative error (%RE). Trueness, a related CLSI term, is the closeness of agreement between the average of a number of replicate measured quantity values and a reference quantity value.
- b. **Analyte:** Component represented in the name of a measurable quantity. The closely related term measurand is defined as the particular quantity subject to measurement.
- c. **Analytical sensitivity:** There are several alternative uses of this term. Most commonly, and for the purposes of this Validation Plan, it is used interchangeably with limit of detection. It is also used to describe the ability of an analytical method to assess small variations of the concentration of an analyte, such as the slope of the calibration curve (IUPAC).
- d. **Analytical specificity:** Ability of a test or procedure to correctly identify or quantify an entity, including in the presence of interfering substance(s) or phenomena.
- e. **Calibration:** Set of operations that establish, under specified conditions, the relationship between values of quantities indicated by a measuring instrument or measuring system, or values represented by a material measure or a reference material, and the corresponding values realized by standards. Under CLIA, calibration refers to the process of testing and adjusting an instrument, kit, or test system, to provide a known relationship between the measurement response and the value of the substance being measured by the test procedure (42 CFR 493.1217).
- f. **Calibrator:** A substance, material, or article intended to be used to establish the measurement relationships of a diagnostic medical device.
- g. **CLIA:** Clinical Laboratory Improvement Amendments of 1988. Congressional legislation that defined and requires specific quality assurance practices in clinical laboratories.
- h. **CLSI:** Clinical and Laboratory Standards Institute.

**Theranos Confidential**

Page 4 of 15

Any retransmissions, reproductions, dissemination or other use of these materials by persons or entities other than the intended recipient is prohibited. This document supersedes all earlier or previous documents unless approved in writing.

TMP-00009 Rev. A, Released 08/01/13

<b>theranos</b>	LDT Validation Report	Theranos Total Protein Assay	Rev:
		<b>CL-RPT-14038</b>	<b>1</b>
Description	Validation Report for Modified Siemens Assay of Total Protein in Lithium Heparin Plasma		
Originator: Curtis Schneider		Date: 09/24/2013	

- i. **Coefficient of Variation:** The ratio of the standard deviation to the average, often multiplied by 100 and expressed as a percentage, abbreviated as %CV .
- j. **Colorimetry:** A technique used to determine the concentration of colored compound(s) in solution.
- k. **Interfering substance:** A substance or quantity thereof that is not the measurand but that affects the result of the measurement.
- l. **IUPAC:** International Union of Pure and Applied Chemistry
- m. **LDT:** Laboratory –developed Test.
- n. **Linearity:** Linearity is the ability of a quantitative analytical method to provide results that are directly proportional to the concentrations of an analyte in test samples, within a given measuring interval. It is an important parameter to confirm when evaluating an analytical method because it verifies correct interpolation of results between points.
- o. **LMR:** Lower end of the measuring range is the lowest level at which defined conditions, including all stated characteristic of the method, are met.
- p. **LoB:** Limit of Blank is the highest value in a series of results on a sample that contains no analyte
- q. **LoD:** Limit of Detection is the lowest amount of analyte in a sample that can be detected with stated probability, although perhaps not quantified as an exact value.
- r. **LoQ:** When used without a prefix, the Limit of Quantitation is the lowest actual concentration at which an analyte is reliably detected and at which uncertainty of the test result is less than or equal to the goal set by the manufacturer or laboratory. The term may also be used with prefixes L for lower (LLOQ) and U for upper (ULOQ), respectively. Note:  $LoB < LoD \leq LoQ$ .
- s. **Matrix:** All components of a material system, except the analyte. A specimen matrix is the biological milieu in which an analyte exists (e.g., plasma, serum, urine, or other body fluids).

**Theranos Confidential**

**Page 5 of 15**

Any retransmissions, reproductions, dissemination or other use of these materials by persons or entities other than the intended recipient is prohibited. This document supersedes all earlier or previous documents unless approved in writing.

TMP-00009 Rev. A, Released 08/01/13

<b>theranos</b>	LDT Validation Report	Theranos Total Protein Assay	Rev:
		<b>CL-RPT-14038</b>	<b>1</b>
Description	Validation Report for Modified Siemens Assay of Total Protein in Lithium Heparin Plasma		
Originator: Curtis Schneider		Date: 09/24/2013	

- t. **Measuring Interval (reportable range; analytical measurement range or AMR):**  
A measuring interval consists of all numeric values between the lower and upper numeric values for which a method can produce quantitative results suitable for clinical use. Where applicable, a linearity study is frequently used to establish or verify the measuring interval that can be reported for a measurement method. Alternatively, the lower limit of the measuring interval may be assigned as the Lower Limit of Quantification (LLOQ).
- u. **Precision:** Precision is the closeness of agreement between indications or measured quantity values obtained by replicate measurements on the same or similar objects under specified conditions. It is usually expressed numerically in terms of standard deviation (SD) or percent Coefficient of Variation (%CV).
- v. **Reference interval:** The interval between and including two reference limits. It is common practice to define a reference limit so a stated fraction of the reference values is less than or equal, or greater than or equal, to the respective upper or lower limit.
- w. **SOP:** Standard Operating Procedure.
- x. **Spectrophotometry:** The quantitative measurement of the transmission (or reflection) properties of a material as a function of wavelength.
- y. **Testing System:** The entirety of the testing process, including instrument, sample, reagents, supplies, and procedures. Personnel are sometimes included in the definition.

**Theranos Confidential**

**Page 6 of 15**

Any retransmissions, reproductions, dissemination or other use of these materials by persons or entities other than the intended recipient is prohibited. This document supersedes all earlier or previous documents unless approved in writing.

TMP-00009 Rev. A, Released 08/01/13

<b>theranos</b>	LDT Validation Report	Theranos Total Protein Assay	Rev:
		<b>CL-RPT-14038</b>	<b>1</b>
Description	Validation Report for Modified Siemens Assay of Total Protein in Lithium Heparin Plasma		
Originator: Curtis Schneider		Date: 09/24/2013	

#### IV. Pre-clinical Validation

##### a. Analytical Measurement Range

###### i. Limits of Blank, Detection and Quantitation

The limits of blank, detection, and quantitation determined were 0.00 g/dL, 0.09 g/dL, and 0.96 g/dL (96% recovery), respectively.

###### Limit of blank

CLSI guideline EP17-A section 4.3.1

Level	Number of samples	N	Mean	SD
Blank	1	20	0.00	0.00
Alpha	5%			
Parametric LoB	0.00			

###### Limit of detection

CLSI guideline EP17-A section 4.3.2

Level	Number of samples	N	Pooled SD
Low	1	20	0.05
Beta	5%		
Parametric LoD	0.09		

###### Limit of quantitation

CLSI guideline EP17-A section 5.1

Level	Number of samples	N
Low	1	20
Bias	-0.04	
Pooled imprecision	0.05	
95% total error	-0.15	
Allowable error	-	
LoQ	0.09	

**Theranos Confidential**

Page 7 of 15

Any retransmissions, reproductions, dissemination or other use of these materials by persons or entities other than the intended recipient is prohibited. This document supersedes all earlier or previous documents unless approved in writing.

TMP-00009 Rev. A, Released 08/01/13

<b>theranos</b>	LDT Validation Report	Theranos Total Protein Assay	Rev:
		<b>CL-RPT-14038</b>	<b>1</b>
Description	Validation Report for Modified Siemens Assay of Total Protein in Lithium Heparin Plasma		
Originator: Curtis Schneider		Date: 09/24/2013	

Level	Sample	n	Assigned value	Mean	Median	SD	CV
Blank	1	20	0	0.00	0.00	0.00	-
Low	1	20	1	0.96	1.00	0.05	5.3%

### ii. Linearity

The Analytical Measurement Range (AMR) including linear measurement interval has been determined by Siemens. Refer to the **Analytical Range** section of the manufacturer product information insert for additional details.

### b. Analytical Specificity

The analytical specificity for this assay was determined by observing the effect of triglycerides (400 mg/dL) on the recovery of total protein (6.92 g/dL) in a spiked plasma sample. No significant interference (NSI) was determined if the mean analyte concentration of an interferent-spiked sample reported within 10% of the mean analyte concentration of an un-spiked sample. Recovery of total protein in the presence of triglycerides was 99.4% (see table below).

Analyte:	Interferent:	
Total Protein (g/dL)	Triglycerides (mg/dL)	% Recovery
6.92	400	99.4*

\* NSI observed at interferent level tested.

### c. Precision

**Theranos Confidential**

Page 8 of 15

Any retransmissions, reproductions, dissemination or other use of these materials by persons or entities other than the intended recipient is prohibited. This document supersedes all earlier or previous documents unless approved in writing.

TMP-00009 Rev. A, Released 08/01/13



<b>theranos</b>	LDT Validation Report	Theranos Total Protein Assay	Rev:
		<b>CL-RPT-14038</b>	<b>1</b>
Description		Validation Report for Modified Siemens Assay of Total Protein in Lithium Heparin Plasma	
Originator: Curtis Schneider		Date: 09/24/2013	

**Level = L1**

Number of observations 80  
 Number of runs 40  
 Number of days 20  
 Runs per day 2  
 Replicates per run 2

Mean 4.04

	SD	95% CI	CV
Repeatability	0.08	0.07 to 0.11	2.1%
Between-run	0.00		0.0%
Between-day	0.09		2.1%
Within-laboratory	0.12	0.10 to 0.15	2.9%

**Level = L2**

Number of observations 80  
 Number of runs 40  
 Number of days 20  
 Runs per day 2  
 Replicates per run 2

Mean 5.84

	SD	95% CI	CV
Repeatability	0.07	0.06 to 0.09	1.2%
Between-run	0.03		0.5%
Between-day	0.07		1.2%
Within-laboratory	0.10	0.08 to 0.13	1.7%

**Theranos Confidential**

**Page 9 of 15**

Any retransmissions, reproductions, dissemination or other use of these materials by persons or entities other than the intended recipient is prohibited. This document supersedes all earlier or previous documents unless approved in writing.

TMP-00009 Rev. A, Released 08/01/13

<b>theranos</b>	LDT Validation Report	Theranos Total Protein Assay <b>CL-RPT-14038</b>	Rev: <b>1</b>
	Validation Report for Modified Siemens Assay of Total Protein in Lithium Heparin Plasma		
Description		Validation Report for Modified Siemens Assay of Total Protein in Lithium Heparin Plasma	
Originator: Curtis Schneider		Date: 09/24/2013	

Level = L3

Number of observations	78
Number of runs	39
Number of runs excluded	1
Number of days	20
% of days with 1 run	5%
Runs per day	2
Replicates per run	2

CLSI guideline EP05-A2 section 10.4 recommends a minimum of 40 runs

Mean	7.26		
	SD	95% CI	CV
Repeatability	0.11	0.09 to 0.14	1.5%
Between-run	0.00		0.0%
Between-day	0.05		0.7%
Within-laboratory	0.12	0.10 to 0.14	1.6%

## V. Clinical Validation

### a. Method Comparison with Predicate (Accuracy/Comparability)

To test the accuracy of the assay on the Theranos System, thirty nine (39) unique patient samples were screened on the predicate method (Siemens, Advia) and on the Theranos method. Using the predicate method twenty (26) values were within the reference range (5.7 – 8.2 mg/dL), six (6) were below the reference range, and seven (7) were above the reference range. Based on the results of the data examination, either a simple linear regression or alternative procedures were used to estimate expected (average) bias and the confidence interval of expected bias at the desired medical decision level(s) as per CLSI guidance EP09-A2. StatisPro was used for bias calculations. These estimates were compared with internal criteria to judge the acceptability of the Theranos method. Each sample was run in duplicate on the predicate, and the average used for comparison to the Theranos method. Some samples were stored before analysis on both methods. If the confidence interval for the predicted bias includes the defined acceptable bias or if the acceptable bias is greater than the higher limit of the confidence interval of the predicted bias, then the data do not show that the bias of the Theranos method is different from the acceptable

**Theranos Confidential**

Page 10 of 15

Any retransmissions, reproductions, dissemination or other use of these materials by persons or entities other than the intended recipient is prohibited. This document supersedes all earlier or previous documents unless approved in writing.


TMP-00009 Rev. A, Released 08/01/13

<b>theranos</b>	LDT Validation Report	Theranos Total Protein Assay	Rev:
		<b>CL-RPT-14038</b>	<b>1</b>
Description	Validation Report for Modified Siemens Assay of Total Protein in Lithium Heparin Plasma		
Originator: Curtis Schneider		Date: 09/24/2013	

bias or there is a high probability (97%) that the predicated bias is acceptable, respectively. The acceptable bias at each medical decision level was determined based on the total allowable error (TEa) minus the measured precision at the level closest to that decision level. Total allowable error (TEa) was taken from American Proficiency Institute (API) peer proficiency testing criteria or CLIA proficiency testing criteria for acceptable analytical performance, as printed in the Federal Register February 28, 1992;57(40):7002-186, when available. The TEa for Total Protein is 10%. The table below shows the allowable bias and precision at 2 levels (values shown in parentheses) and the corresponding closest medical decision limits.

Table 1. Allowable Bias and Precision at the Medical Decision Levels

<b>Medical Decision Levels (g/dL)</b>	<b>4.5 (4.0)</b>	<b>6.0 (5.8)</b>
<b>Precision (%)</b>	<b>2.1</b>	<b>1.2</b>
<b>Allowable Bias (%)</b>	<b>7.9</b>	<b>8.8</b>



**Figure 1.** Graph showing Theranos method versus Predicate Method (Siemens Advia). Simple linear regression was used to establish a slope, intercept and an r<sup>2</sup>. The slope, intercept and clinical correlation were determined to be 1.04, -0.18 and 0.99 respectively.

Any retransmissions, reproductions, dissemination or other use of these materials by persons or entities other than the intended recipient is prohibited. This document supersedes all earlier or previous documents unless approved in writing.

<b>theranos</b>	LDT Validation Report	Theranos Total Protein Assay	Rev:
		<b>CL-RPT-14038</b>	<b>1</b>
Description	Validation Report for Modified Siemens Assay of Total Protein in Lithium Heparin Plasma		
Originator: Curtis Schneider		Date: 09/24/2013	

**Comparability**

CLSI guideline EP09-A2-IR section 7


Level ID	Value	Difference	SE	95% CI	Allowable difference
	4.5	-0.13	0.070	-0.27 to 0.01	0.32
	6.0	-0.10	0.055	-0.21 to 0.01	0.43
	8.0	-0.06	0.053	-0.17 to 0.04	0.57

Difference is less than allowable bias: 7.15%.

The difference between the two methods is not greater than the allowable difference. The performance requirement is verified.

**b. Transference and Verification of Reference Interval (Venous)**

Reference ranges were modified by applying the regression equation to the lower and upper reference limits of existing reference interval to generate a new reference range. New reference ranges were verified with venous samples using forty (40) new normal subjects. For a reference range to pass verification, 95% of values should fall within the upper and lower reference limits and 5% or fewer values fall outside of the upper and lower reference limits. For venous verification 38 (95%) values fell within the new reference range and 2 (5%) values fell outside the new reference range. See graph below for venous samples verification.




**Figure 2.** Graph showing venous sample reference range verification.

<b>theranos</b>	LDT Validation Report	Theranos Total Protein Assay	Rev:
		<b>CL-RPT-14038</b>	<b>1</b>
Description	Validation Report for Modified Siemens Assay of Total Protein in Lithium Heparin Plasma		
Originator: Curtis Schneider		Date: 09/24/2013	

**c. Verification of Reference Interval with Finger Stick Samples**

New reference ranges were also verified with venous matched finger sticks from thirty five (35) new normal subjects. For finger stick verification 34 (97.1%) values fell within the new reference range and 1 (2.9%) values fell outside the new reference range. See graphs below for finger stick samples verification.



**Figure 3.** Graph showing Finger stick sample reference range verification.

Theranos.Analyte	Anti-coagulant	Existing Reference Range (g/dL)		New Reference Range (g/dL)	
		Reference Range (low)	Reference Range (High)	Transferred RR (low)	Transferred RR (high)
Total Protein	Lithium Heparin	5.7	8.2	5.7	8.3

The new reference range for finger stick Total Protein was determined to be 5.7 – 8.3 g/dL.

<b>theranos</b>	LDT Validation Report	Theranos Total Protein Assay	Rev:
		<b>CL-RPT-14038</b>	<b>1</b>
Description	Validation Report for Modified Siemens Assay of Total Protein in Lithium Heparin Plasma		
Originator: Curtis Schneider		Date: 09/24/2013	

## VI. Stability

### a. Reagents

#### On-board Reagent Stability

System	Stability without Reagent Container Inserts*	Stability with Reagent Container Inserts
ADVIA 1200	21 days	60 days
ADVIA 1650/1800	30 days	60 days
ADVIA 2400	30 days	60 days

\*If your test workload requires the longer on-board stability, use a reagent container inserts (REF 02991886) in R1 and R2 containers.

For all systems, unopened reagents are stable until the expiration date printed on the product label when stored at 15° - 25°C. Do not freeze reagents.

For complete details, refer to the Methods Introduction section of the system-specific Operator's Guide.

### b. Sample

Plasma samples for Total Protein analysis are stable for 2 weeks at 2-8 °C, or at least 90 days at -20 °C.

### c. Calibrators

Siemens Chemistry Calibrators should be stored at 2-8 °C, protected from light, and are stable until the expiration date on the vial label. Opened calibrators are stable for 48 hours, except for total and direct, which are stable for 8 hours.

## REVISION HISTORY

Revision Level	Effective Date	Initiator	ECO Number
A	11/06/2013	A. Rosendorff	CL ECO-00117

**Theranos Confidential**

Page 14 of 15

Any retransmissions, reproductions, dissemination or other use of these materials by persons or entities other than the intended recipient is prohibited. This document supersedes all earlier or previous documents unless approved in writing.

TMP-00009 Rev. A, Released 08/01/13

<b>theranos</b>	LDT Validation Report	Theranos Total Protein Assay	Rev:
		<b>CL-RPT-14038</b>	<b>1</b>
Description	Validation Report for Modified Siemens Assay of Total Protein in Lithium Heparin Plasma		
Originator: Curtis Schneider		Date: 09/24/2013	

Section Number	Description and Justification of Changes
All	Initial Release

Theranos Internal Only

Any retransmissions, reproductions, dissemination or other use of these materials by persons or entities other than the intended recipient is prohibited. This document supersedes all earlier or previous documents unless approved in writing.

TMP-00009 Rev. A, Released 08/01/13