

Surekha Gangakhedkar

From: Daniel Young
Sent: Friday, August 30, 2013 7:13 PM
To: Surekha Gangakhedkar; Elizabeth Holmes
Cc: Sharada Sivaraman
Subject: RE: missing data from 3.5 runs

I recommend repeating the failed runs. We will note in the study report that x number of failed runs (due to the wrong cartridges being used) were rerun.

Daniel

From: Surekha Gangakhedkar
Sent: Friday, August 30, 2013 1:39 PM
To: Elizabeth Holmes; Daniel Young
Cc: Sharada Sivaraman
Subject: FW: missing data from 3.5 runs

Hi Elizabeth,

We had issues with all runs on the 3.5 yesterday as none of the runs completed due to temperature not being reached and pipette initialization errors so we have to repeat those runs (for Vitamin D). The software team is working on addressing other issues.

Daniel – Please review data so far on the 3.0 for TPSA. We did have lots of errors due the cartridges failing due to the blood tip drop off issue (due to the absence of labels). We are presently testing 30 cartridges with a fix in place. Please advise on if the data points will have to be repeated ?

Thanks,
Surekha

From: Sekhar Variam
Sent: Friday, August 30, 2013 1:25 PM
To: Tina Noyes; Earl Solis; Samartha Anekal; Michael Craig
Cc: Surekha Gangakhedkar; Sharada Sivaraman
Subject: RE: missing data from 3.5 runs

Hi Tina:

None of the runs completed and we don't have any results on the database.

Protocol Execution failed for the following barcodes because temperature was not reached

91245647700100442
91245647700100498
91245647700100500

For the remaining barcodes, I think there was an issue in initializing the pipette. I confirmed it on E000240, but could not confirm on the remaining devices.