

Message

From: Sunny Balwani [/O=Theranos Organization/OU=First Administrative Group/CN=Recipients/CN=SBALWANI]
Sent: 8/29/2013 5:14:20 PM
To: Elizabeth Holmes [eholmes@theranos.com]
Subject: RE: Concerns about the launch

We will talk to him.

I assume we can use Siemens chemistry for glucose and sodium? If yes then we should

As for operational, we should tell him we will limit # of samples to 50 per day for 2-4 weeks until we have addressed all staffing and training needs.

From: Elizabeth Holmes
Sent: Thursday, August 29, 2013 10:11 AM
To: Sunny Balwani
Subject: Fwd: Concerns about the launch

Begin forwarded message:

From: Adam Rosendorff <arosendorff@theranos.com>
Date: August 29, 2013, 10:06:10 AM PDT
To: Elizabeth Holmes <eholmes@theranos.com>
Cc: Daniel Young <dyoung@theranos.com>
Subject: Concerns about the launch

Hi Elizabeth

I have some medical and operational concerns about our readiness for 9/9.

1) Medical

GLUCOSE

Significant negative bias exists for glucose, 95% confidence interval for bias (negative):

Medical Decision Level: 45 mg/dL 120 mg/dL 180 mg/dL

11%

-20% to -12% -16% to -11% -16% to -

N=123

This is again reflected in the reference range which is coming out as 20-97 mg/dL on the Theranos chemistries whereas on the Siemens chemistries we are establishing a de-novo reference range of 56-117 mg/dL, which compares nicely with the published range of 70-110 which means we are doing our studies properly. The recommended maximum bias is 5% Our CV% is also a little high (8%). The CLIA total allowable error (precision plus bias) is 10%, which means that our entire error budget has already been spent just on precision. I am not sure how comfortable I am applying a correction factor for glucose because criteria for diabetes for instance, have been set at a fasting blood glucose of 126 mg/dL, which is fairly standard across analyzers without post-hoc corrections. Maybe we need to work on our calibrators or chemistry?

SODIUM

I noticed that for one of the demos a few patients had sodium values in the 120s, whereas 99% of healthy patients should have sodiums between 135-145 mEq/L. The allowable bias for sodium should be 2-3 mEq/L

2) Operational

We are understaffed and not trained on the new protocols. We would need at least 2 CLSes downstairs in Normandy to do the pipetting from the BCDs/ handling of the TECAN, and 2 operators upstairs to do the QC and ADVIA operation, in addition to a "float" CLS that can fill in as needed. This means we need 5 CLSes and we only have 3. I am working with Curtis and Chinmay to finalize SOPs, and with Erez to crunch through the validation data.

I would like us to be the best that we can be. A few more weeks to sort through these medical and logistical issues, and getting the proper level of training and staffing would help us tremendously.

Regards,

Adam

Adam Rosendorff, MD, FASCP

Theranos, Inc

(650) 856-4412 (Office)

(650) 823-4953 (Mobile)

(650) 852-9594 (Fax)

