

HIV-2 Assay Development Report

Theranos, Inc.

This Development Report contains Theranos Confidential Information. Any further dissemination, use or disclosure of the Report, in whole or in part, is strictly prohibited.

Table of Contents

[TOC at "Heading 1,1]

Theranos Internal Only

[TOC \h \z \c "Table"]

Theranos Internal Only

1 Purpose

This document describes the studies conducted during the development of the HIV-2 IgG assay.

2 Scope

This report includes relevant experiments in the development of the HIV-2 IgG assay.

3 Definitions

Term	Definition
HIV-2	Human Immunodeficiency Virus Type 2
RF	Rheumatoid factor
RLU	Relative light unit
IgG	Immunoglobulin G
HAMA	Human anti-mouse antibodies
NPA	Negative Percent Agreement
PPA	Positive Percent Agreement

Table [SEQ Table * ARABIC]: Definitions

4 Assay principle

The Theranos HIV-2 assay is an in vitro diagnostic immunoassay for the qualitative determination of IgG antibodies to HIV-2 in human serum, plasma and whole blood. In this assay, the sample is incubated with a solid surface coated with HIV-2 antigens. Antigen-antibody complexes will form if the patient sample contains anti-HIV-2 antibodies. After a number of washing steps to remove unbound sample, the solid surface is incubated with an HIV-2 peptide conjugated to a reporter enzyme. After a number of washing steps to remove unbound detection conjugate, the solid surface is incubated with a chemiluminescent substrate and RLUs are measured using the Theranos Analyzer. The RLU data is used to calculate each sample's index value by comparison with on-board calibrators.

5 Reference Assay

- 5.1 The following assay was used as reference method:
HIV-2, ADVIA Centaur XP

6 Preliminary clinical study

- 6.1 The goal of this experiment is to evaluate the Theranos assay for its ability to detect anti-HIV-2 antibody in a group of individuals.
- 6.2 Procedure: A total of 281 patient samples, including individuals at high risk, pregnant individuals and individuals with signs and symptoms of HIV-2 infection, were analyzed using the Theranos system and compared with results from the reference method, ADVIA Centaur XP HIV-2 IgG assay.
- 6.3 Acceptance criteria: Positive percent agreement (sensitivity) should be 100% and negative percent agreement (specificity) should be 100%.
- 6.4 Results: The overall positive percent agreement was 100% and negative percent agreement was 100%.

		High Risk, including pregnancy					
		Theranos			Total		
Centaur HIV-2 IgG		Positive	Negative	Equivocal	Total	100.00%	NPA
	Positive	0	0	0	0		
	Negative	0	119	0	119		
	Equivocal	0	0	0	0		

		Symptomatic subjects + Positive samples					
		Theranos			Total		
Centaur HIV-2 IgG		Positive	Negative	Equivocal	Total	100.00%	PPA
	Positive	36	0	0	36		
	Negative	0	0	0	0		
	Equivocal	0	0	0	0		

		Additional panels					
		Theranos			Total		
Centaur HIV-2 IgG		Positive	Negative	Equivocal	Total	100.00%	PPA
	Positive	15	0	0	15		
	Negative	0	111	0	111		

	Equivocal	0	0	0	0
--	-----------	---	---	---	---

All samples

		Theranos					
Centaur HIV-2 IgG		Positive	Negative	Equivocal	Total		
	Positive	51	0	0	51	100.00%	PPA
	Negative	0	230	0	230	100.00%	NPA
	Equivocal	0	0	0	0		

*Table [SEQ Table * ARABIC]: Clinical Study Results*

6.5 Conclusions: Our results meet the acceptance criteria and additional samples (10000 samples) will be run during validation, which will further indicate and support that the assay meets the acceptance criteria

7 Analytical specificity (cross-reactivity)

- 7.1 The goal of this experiment is to test the Theranos HIV-2 assay for cross reactivity to medical conditions unrelated to HIV-2.
- 7.2 Procedure: 201 samples from 17 cross-reactive groups were tested using the HIV-2 Theranos Analyzer and compared to ADVIA Centaur.
- 7.3 Acceptance criteria: minimal cross reactivity to medical conditions unrelated to HIV-2. Positive % agreement should be 100% and negative % agreement should be 100%
- 7.4 Results: 100% of negative samples were found to be non-reactive to both the Theranos and Reference method (negative % agreement was 100%), and 100% of positive samples were found to be reactive for both assays (Positive % agreement was 100%). Summary is shown in the following table:

Category #	Cross Reactant	Samples Tested	In Agreement
1	Elevated IgG	14	14
2	ANA	14	14
3	EBV IgG	14	14
4	RF+	14	14
5	CMV	14	14
6	Toxoplasma Gondii	14	14
7	HBV	14	14
8	HBsAg	7	7
9	Hemodialysis	12	12

10	Graves Disease	12	12
11	Flu	10	10
12	Scleroderma	14	14
13	Rubella	10	10
14	Lupus	10	10
15	HSV 1/2	12	12
16	Elevated IgM	10	10
17	HCV	6	6

Table [SEQ Table * ARABIC], cross-reactivity to different medical conditions

7.5 Conclusions: No cross-reactivity observed. Acceptance criteria for cross-reactivity was met

8 Interference

8.1 The goal of this experiment was to test the effect of lipemic, hemolyzed, icteric, high-protein, RF+, and HAMA+ samples on Theranos assay

Procedure: An HIV-2 positive and an HIV-2 negative sample were spiked with 750 mg/dL triglycerides, 1000 mg/dL hemoglobin, 20 mg/dL bilirubin, 2g/dL total protein, 500 ng/mL HAMA+, or 2000 IU/mL RF+. These spiked samples were tested on the Theranos Analyzer and compared with non-spiked positive and negative HIV-2 samples to determine percent interference. Percent interference was calculated with the equation below. x_M is the back calculated measured value of the spiked sample on the Theranos assay and x_T is the value of the unspiked sample on the Theranos assay.

- If reference COI < 0.5, absolute value of the interference must be less than 0.1 COI, i.e.

$$\text{Interference} = x_M - x_T$$

- If reference COI > 0.5, absolute value of the interference must be less than 20%, i.e.

$$\text{Interference} = (x_M - x_T) / x_T, \%$$

- 8.2 Acceptance criteria: The absolute value of the percent interference must be less than 20%.
- 8.3 Results: All spiked samples showed no difference in diagnosis between Sample and Control. Summary is shown in the following tables:

Sample Type	Interfering Substance	Interference
HIV-2 Positive (COI ~ 1.5)	Hemoglobin (1000 mg/dL)	2%
	Bilirubin (20 mg/dL)	-14%
	Triglycerides (750 mg/dL)	1%
	Total Protein (2 g/dL)	10%
	RF+ (2000 IU/mL)	14%
	HAMA+ (500 ng/mL)	7%
HIV-2 Negative (COI = 0.02)	Hemoglobin (1000 mg/dL)	0.01
	Bilirubin (20 mg/dL)	0.01
	Triglycerides (750 mg/dL)	0.02
	Total Protein (2 g/dL)	0.01
	RF+ (2000 IU/mL)	0.01
	HAMA+ (500 ng/mL)	0.00

Table [SEQ Table 1*ARABIC]: Interference to endogenous compounds

8.4 Conclusions: No interference was observed. The acceptance criteria was met.

9 Matrix comparison

- 9.1 The goal of this experiment was to use Theranos assay to compare Venous EDTA Plasma with Venous Lithium Heparin Plasma and Venous Serum.
- 9.2 Procedure: 14 unique patients each donated 1 Serum tube, 1 K₂EDTA tube, and 1 Li-Hep tube. Plasma from all BCDs from the same patient were pooled. Due to the fact that all in-house collected samples were HIV-2 negative, the matched samples from 14 patients were spiked using different amounts of an HIV-2 positive serum sample. HIV-2 levels were determined by the Theranos method from all 3 matrices and compared to each other.
- 9.3 Acceptance criteria: Each EDTA or Lithium Heparin sample must be within +/- 30% difference from the serum sample.
- 9.4 Results: Table 5 summarizes the matrix comparison results for the HIV-2 assay. One Lithium-Heparin sample is discrepant when compared to the serum reference.

Sample	Venous EDTA	Venous Li-Hep
1	12.4	-6.9
2	-14.6	-5.8
3	-26.2	-5.2
4	-21.6	-8.0
5	-3.9	-11.9
6	-23.1	-20.1
7	-5.1	-1.2
8	-7.3	2.6
9	-27.8	-22.9
10	-3.1	5.3
11	-12.4	-8.3
12	-25.5	18.3
13	-28.9	31.3
14	-21.8	-11.7

Table [SEQ Table * ARABIC]: Matrix Comparison Results

9.5 Conclusions: Only one discrepant sample compared with the venous serum reference demonstrating matrix equivalence for the HIV-2 IgG assay.

10 Conclusions

We successfully developed an ELISA sandwich assay to detect HIV-2 IgG in different matrices. The assay meets all requirements for sensitivity, specificity, cross reactivity, interference, and matrix effect.

Theranos Internal Only